Volume 29, Issue 2, March 2023, Page (620-625)



DOI

https://dx.doi.org/10.21608/zumj.2022.161051.2634

Manuscript ID ZUMJ-2209-2634 (R1)

10.21608/zumj.2022.161051.2634

ORIGINAL ARTICLE

Serum Chitinase3-like Protein 1 (YKL-40) in Wheezy Infants. El Sayed Wagdy¹; Ibrahim Dina Ahmed²

1Pediatrics Department, Al Ahrar Teaching Hospital, Zagazig, Egypt 2Clinical Pathology Department, Al Ahrar Teaching Hospital, Zagazig, Egypt

Corresponding author El Sayed, Wagdy E-mail: wagdyattia@yahoo.com

 Submit Date
 2022-09-06

 Revise Date
 2022-10-10

 Accept Date
 2022-11-02

ABSTRACT

Background: Chitinase-3-Like Protein 1(YKL40) is an inflammatory biomarker secreted from neutrophils, epithelial cells and smooth muscle cells of airways. It has been expressed in diseases characterized by inflammation and tissue remolding of respiratory tract e.g., asthma. Little reports were described in wheezing infants. The aim of this study is to estimate serum YKL40 in wheezy infant and study its role in pathogenesis in different types of wheezes which may help in diagnosis and management of the disease.

Methods: The study was carried on 49 infants with wheezes ,23 has recurrent wheezes ,26 has no recurrent wheezes. Detailed medical history and examination were done. Laboratory investigations were done including Complete blood count (CBC) with differential counts, C-reactive protein (CRP), blood gases were done. Chest x-ray and C.T. scan were done when required for persistent wheezing, YKL-140 was assessed by enzyme linked immune assay kits.

Results: Serum YKL 40 did not show any relation to neonatal intensive care unit (NICU) admission or recurrence of wheezing. Cases with family history of asthma showed higher level compared with those without family history. Also, critical cases showed higher values compared with stable cases. However, cases without family history of wheezing and stable cases showed no difference in serum YKL 40 compared to control. No significant difference between cases and control regarding WBCs, RBCs, platelets and CRP, while cases showed lower Hb and higher, neutrophil, serum YKL40 compared with control. Serum YKL140 showed positive correlation with neutrophil counts.

Conclusion: Serum YKL140 was significantly higher in wheezing infants compared to control. Cases with family history of asthma showed higher level of this marker compared with those without family history of asthma. It showed positive correlation with blood neutrophil. These



results have diagnostic and prognostic significance regarding infant wheezes.

Keywords: Serum Chitinase 3-like protien1 KYL-40, infant wheezes, atopic, recurrent, non-recurrent.

INTRODUCTION

S erum chitinase-3 like protien-1 "YKL-40" is an inflammatory biomarker associated with disease activity and mortality in patients with diseases characterized by inflammation and tissue remodeling [1]. The site of formation, action and relation to disease state of YKL-40 is subject of research in recent years [2]. It is secreted from cells in the airway mucosa, i.e., neutrophils, macrophages, and airway epithelial cells [3]. It is induced at sites of inflammation, infection, and tissue remodeling, suggesting that this protein plays active roles in anti-infective defense and repair responses [4].

In adults, increased serum YKL-40 has been associated with asthma, and its levels correlate with increasing disease severity and bronchial wall thickness as well as decreasing lung function [5].

In pediatrics, YKL-40 it was studied in different pulmonary diseases for monitoring and diagnosis. YKL-40 in bronco-alveolar lavage is reflective of inflammation in early cystic fibrosis lung disease in children [6]. Significant association between serum chitinase-3-like protein-1 (YKL-40) and lung function in children with asthmatic symptoms was detected.

This association was focused on indices that reflect small airway function [7]. It was found to be higher in children with pneumonia and has prognostic value in those with viral pneumonia [8].

Wheezing is a common challenge in early childhood. About 50 % of children experienced one wheezing episode in 1st 6 ys. One third had one episode of wheezing in the preceding 6 winter months [9]. There are 3 phenotypes of wheezes: 1-Early transient wheeze (20%) start before 3 years and ends before 6 ys 2-Late onset (15%) starts between 3-6 ys , and 3-Persistent wheezes (14%) starts before 3 ys and continue after 6 years and divided into 2 types atopic or IgE mediated and non-atopic type [10].

Little reports available about serum YKL 40 in infant wheezes. In preschool children YKL 40 was significantly higher in wheezing infant compared to control but not related to persistent wheezing [11]. It was higher in post infectious bronchiolitis obliterans exacerbation compared to acute bronchiolitis [12].

For further information about its role of in pathogenesis and management in respiratory diseases in children we measured YKL-40 in wheezy infants. Also, we hope to study its role in different types of wheezes.

METHODS

Patients group: Forty-nine patients who suffered from wheezing were included in this study. 23 were suffered from recurrent wheezes while 26 has no recurrent wheezing.

Detailed medical history including admission in NICU, chest diseases, allergies and smokes exposure. Family history of asthma, atopic and immunological diseases were taken. Complete physical examination and growth assessment were done including weight, height. Cases were classified according to clinical presentation into two groups :-1-stable cases those cases presented with feeding amount less than usual but still > half of normal amount, signs of breathing difficulty, respiratory rate (R.R.) <60/min , no grunting , moderate retraction, wheezing, O2 saturation 88-92% and generally look irritable 2-ceitical(severe cases), those cases presented as moderate cases + grunting, severe retraction, O2 saturation <88%, feeding decreased to less than half amount, generally look lethargic [13] Laboratory investigation included CBC with differential count, CRP. Blood gases was assessed. Chest xray and C.T. scan were done when required for persistent wheezing, YKL-40 was assessed by enzyme linked immune assay kits.

Sample Collection: Blood samples were obtained by venipuncture from each subject and divided in two tubes. 1st tube with K2EDTA for blood count and blood film 2nd tube is allowed to coagulate for 30 minutes at room temperature. Serum was separated by centrifugation for 10 min at 3000 rpm and stored at -70 till tests were done.

Quantitation of human chitinase-3 like-proten-1(YKL40): Human serum chitinase-3-like protein 1 (YKL40) was analyzed according to manufacturer instruction by ELISA using human YKL40 kit (Cat No E2063Hu) supplied by Bioassay Technology Laboratory (China).

Principle of the assay: The kit is enzyme linked immunosorbent assay. The plate has been precoated with human (YKL40) antibody. YKL40 present in the sample was added and bind to antibodies coated on the wells and then biotinylated human YKL40 antibody was added and bound to YKL40 in the sample. Then streptavidin- horseradish peroxidase (HRP) was added and bound to biotinylated YKL40 Ab. After incubation unbound streptavidin-HRP was washed away. Substrate solution was then added and color developed in proportion to the amount of YKL40. The reaction is terminated by addition of acidic stop solution and absorbance is measured at 450 nm.

Inclusion criteria infants with wheezes with age ranges from 1-24 months whether recurrent or nonrecurrent, those who are stable (moderate respiratory distress) or critical (has severe distress). Exclusion criteria: Cases with congenital lung lesions, suspected immune deficiency gastroesophageal reflux disease "GERD", possible cardiac, mediastinal lesion was included in the study.

Control group: Fifteen children age and sex matched were selected as control group. These cases were taken from OPD clinic came for routine check or minor illnesses.

STATISTICAL METHODS

Data were analyzed using SPSS 20 computer program (IBM, Endicott, Broome County, New York, United States). Data were expressed as mean \pm SD for categorized variables. Tests of significance "Chi-square and T tests" and correlation study were done when appropriate. The correlation coefficient method was used to correlate different parameters. P value > 0.05 = statistically insignificant, P \leq 0.05* = statistically significant, P < 0.01**= highly significant.

RESULTS

Forty-nine (49) infants suffered from wheezing were included in this study. 23 had recurrent wheezing while 26 had 1st attack of wheezing ,24 has positive (+ve) family history of asthma while

Elsayed, W., Ibrahim, D

25 cases had negative(-ve) family history, 21 had had NICU admission while 28 had no history of NICU admission (Table 1).

Our cases were age, weight and sex matched with control cases (P was 0.446,0.564 and 0.567 respectively) (Table 2).

Serum YKL-40 did not show any relation to NICU admission, recurrence or severity of wheezing, (P was 0.612, 0.334 and 0.513 respectively). Cases with family history of asthma showed higher level compared with those without family history (P was 0.003) (Table 3). However, cases without family history of wheezing showed no difference in serum YKL 140 compared to control (P was 0.562).

No significant difference between cases and control regarding WBCs, RBCs, platelets and CRP (P was,0.727, 0.224,0.479 and 0.323 respectively), while cases showed lower Hb and higher, neutrophil, serum YKL40 (P was 0.031,0.027 and 0.028 respectively) compared with control (Table 4).

Serum YKL 140 did not show any correlation with age, weight, WBCs, RBCs, platelets, CRP or Hb (P was 0.927, 0.429, 0.842, 936, 0.541, and 0.27 respectively) while showed positive correlation with neutrophil count (P was 0.009) (table 5) (Figure 1).

|--|

Symptoms		Number of cases	percentage
Wheezing	Recurrent	23	46.9%
	Non-recurrent	26	53.1%
Family history of asthma	+ve	24	49%
	-ve	25	51%
NICU admission	+ve	21	42.9%
	-ve	28	58.1%
Severity of wheezing	Stable	34	69.4%
	Critical	15	30.6%

NICU: Neonatal Intensive Care Unit

Table 2: Demographic data for cases and control

	Cases (49)	Control (15)	Р	
Age (months) "MW"	2-24(12)	2-24(12)	0.446	
Weight "T-test"	9.19±2.5	9.5±1.8	0.564	
Sex (chi square) M	22	8	0.576	
F	27	7		

MW: Man-Whitney

Table 3: Serum YKL40 relation to clinical data using Man-Whitney Test

Wheezing recurrence	Recurrent (23)	10-50(19)	P:0.334
	Non-recurrent (26)	10-55(17.5)	
Family history of asthma (FHA)	+ve(24)	12-55(26)	P: 0.003
	-ve(25)	12-50(15)	
NICU admission	+ve(21)	10-50(16)	P: 0.612
	-ve(28)	10-55(18)	
Severity of wheezing	Critical (34)	10-55(20)	P:0.513
	Stable (15)	10-50(17.5)	

YKL-40: Serum Chitinase3-like protein1, NICU: Neonatal Intensive Care Unit Table 4: Laboratory data in wheezy infant versus control

		Patients (49)	Control (15)	р
CRP (mg)"MW"		1-38(12)	12-24(12)	0.323
WBCs (×10 ³) "MV	V''	3-20(8.7)	6-13	0.727
Neutrophil (×10 ³)	"MW"	1.5-6(3)	2.3-3(2.7)	0.027
RBCs (10 ⁶) T-test		4.2±0.684	4.43±0.39	0.224
Hb(gm) T-test		11.73±1.27	12.18±1.13	0.031
Platlets×10 ³ T-test		110-620(350)	250-580(350)	0.479
YKL40(ng/ml)	All cases(49"MW")	10-55(19)	(12-24(15)	.042
	-ve FHA (24) "MW"	12-50(15)	12-24(15)	0.562
	+ve FHA (25) "MW"	12-55(25)	12-24(15)	0.002

Elsayed, W., Ibrahim, D

622 | P a g e

YKL-40: Serum Chitinase3-like protein1, NICU: Neonatal Intensive Care Unit, CRP: C-reactive protein, RBCs: Red Cell Counts, Hb: Hemoglobin, FHA: Family history of asthma

Table 5: Serum	YKL40	Correlation	with	Laboratory	Data

Parameter	R	Р
WBCs	0.029	0.842
Neutrophils	0.371	0.009
Platelets	-0.12	0.836
RBCs	-0.084	0.568
Hb	0.161	0.27
CRP	0.089	0.541




Figure 1: Serum YKL-40 relation to Laboratory Parameters

DISCUSSION

Chitinase -3 like protein-1 YKL-40 is an inflammatory biomarker involved in airway remodeling due to its effect on smooth muscle proliferation [14] and its correlation with subepithelial fibrosis and bronchial wall thickness [15,16] Most reports in pediatric patients were described in asthmatic children and obstructive pulmonary diseases [6,7], while little reports were done in wheezy children as primary presentation. The aim of this study to estimate prognostic and possible diagnostic value of YKL 40 in evaluation of severity and type of wheezing which may be reflected on its management plan. No available reports about this marker in infancy period (1 months -2 years. Our cases were divided clinically according to their recurrence, NICU admission, family history of asthma and severity. Mild cases usually treated in OPD and have not been hospitalized and did not require extensive investigation and follow up and were not included in our study. Our cases were age, sex and weight matched with control group (P was 0.15,0.564 and 0.567 respectively).

Wheezy infants showed lower values of Hb (but still in normal range), neutrophil count, serum YKL40 compared to controls (P was: 0.031,0.00 and 0.028 respectively) Neutrophil count showed positive correlation with YKL-40. Unfortunately, no available study about serum YKL-40 in

wheezing in infancy period. James et al [11] found that preschool children (6-44months) with wheezes has significant higher neutrophil count and serum YKL-40 compared to control. Increased neutrophil in wheezing may results from associated viral pro-inflammatory infection and cytokines associate wheezing [17,18, 19]. In that study [11] serum YKL-40 was significantly higher at acute visit and at 3 months follow up but not at 1 year follow up compared to control. Also, it was increased in infant with acute bronchiolitis and bronchiolitis obliterans compared to control [12]. It was increased in children with asthma, and chronic obstructive pulmonary diseases [20-22]. Positive correlation with neutrophil was found in acute bronchiolitis and preschool children with wheezes also found [11,12]. YKL-40 is was an inflammatory marker released from neutrophils and macrophages at time of infection and inflammation in adults [23,24]. This correlation with neutrophil count possibly reflects relation of the marker to severity of inflammation and clinical consequently condition. Allergen challenges in adults induced increase in YKL 40 in BAL despite unknown trigger did not increase its serum level [25-27]. Correlation of this marker to type and severity of inflammation was observed in previous studies [5,12]. In our study serum YKL140 was higher in cases with family history of asthma compared to cases without family history

(P was 0.003 respectively). This results support that rise in this marker can give clue about type of inflammation. This is supported by lack of significant rise in YKL-40 cases with negative family history of asthma in our study (P was 0.562 respectively). No significant difference in serum YKL40 between critical cases and stable cases (p was 0.513). this results possibly due to small sample size in our study. Knihtilä et al found that serum YKL 140 was linked to small airway function in children with asthmatic symptoms as evidenced by correlation with all baseline impulse oscillometry (IOS) indices and higher value in children with significant higher airway hyper (AHR) and bronchodilator reactivity responsiveness (BDR) [7]. In asthmatic children it was higher than control and higher in acute cases compared with stable cases and correlated positively with severity of acute attack. [28,29] In adults with bronchial asthma it correlated inversely with forced expiratory volume (FEV1) [30]. Also, it increased significantly from mild-moderate to severe asthma and was significantly increased in patients with chronic obstructive lung disease compared with severe asthma [31].

Our study reflects that serum YKL-40 increase significantly in wheezy infant compared to control cases with positive family history of asthma and allergic diseases showed higher values. These results reveal prognostic and diagnostic value of the marker in wheezy infant. Previous studies showed elevation of this marker in other allergic conditions as atopic dermatitis. [32]. Relation of this marker to allergies may be explained by association with genetic mutation CHI3L1. Genetic variation in CHI3L1 influenced serum YKL-40 levels and was associated with increased the risk of asthma, bronchial hyper-responsiveness and reduced lung function [30].

Our study has clinical implication. Elevation of serum YKL-40 in wheezy infant may suspect allergic wheezing and possible later development of asthma. This marker may have role in management of wheezing especially of allergic origin. Inhibition of this marker by k284-6111 reduce skin manifestation in atopic dermatitis in animal model through suppression of lactoferrin and consequently inflammatory cytokines [33].

The limitation of this study: Relatively small sample size. Also, no follow up of wheezing cases to correlate the marker to response of treatment and determine type of wheezes. This may need further study with large number of cases with identification of type of wheezing. and relation to this marker.

CONCLUSION AND RECOMMENDATION

Serum YKL-40 was evaluated in wheezy infants. It was significantly higher in wheezing infants compared to control. Cases with family history of asthma showed higher level of this marker compared with those without family history of asthma. It showed positive correlation with blood neutrophil. These results have diagnostic and prognostic significance regarding infant wheezes. List of Abbreviation: chitinase-3-Like Protein 1(YKL-40), complete blood count (CBC), C reactive protein (CRP), gastro-esophageal reflux disease (GERD), Neonatal intensive care unit Respiratory rate (RR). (NICU) impulse oscilometry (IOM), airway hyper reactivity (AHR},bronchodilator responsiveness (BDR).

REFERENCES

- **1-Schultz NA, Johan sen JS.** YKL-40 a protein in the field of translational medicine: a00000 role as a biomarker in cancer patients? Cancers. 2010; 2: 1453-91.
- **2- Lee CG, Hartl D, Lee GR, Koller B, Matsuura H, Da Silva CA, et al.**, Role of breast regression protein 39 (BRP-39)/chitinase 3-like-1 in Th2 and IL-13induced tissue responses and apoptosis.J Exp Med. 2009; 206: 1149-66.
- **3-Park JA, Drazen JM, Tschumperlin DJ.** The chitinase-like protein YKL-40 is secreted by airway epithelial cells at base line and in response to compressive mechanical stress. J Biol Chem. 2010; 285: 29817-25.
- 4- Lee CG, Hartl D, Lee GR, Koller B, Matsuura H, Da Silva CA, Sohn MH, Cohn L et al: Role of breast regression protein 39 (BRP-39)/chitinase 3like-1 in Th2 and IL-13-induced tissue responses and apoptosis.J Exp Med 2009, 206(5):1149–1166.
- 5- Ober C, Tan Z, Sun Y, et al., Effect of variation in CHI3L1 on serum YKL-40 level, risk of asthma, and lung function. N Engl J Med.2008; 358:1682-1691.
- **6- Fantino E, Gangell CL, Hartl D, Sly** PD Airway, but not serum or urinary, levels of YKL-40 reflect inflammation in early cystic fibrosis lung Disease Pulmonary Medicine 2014, 14⊗28),1-6
- 7- Knihtilä H, Kotaniemi-Syrjänen A, Pelkonen AS ,Savinko T and Malmberg LP et al, Serum chitinase-like protein YKL-40 is linked to small airway function in children with asthmatic symptoms 2019;30(8):803-809
- 8- Yang x and Sheng G, YKL-40 levels are associated with disease severity and prognosis of viral pneumonia, but not available in bacterial pneumonia in children BMC Pediatrics 2018, 18:381,1-8
- **9- Bisgaard H and. Szefler S.,** Prevalence of asthma like symptom in young children Pediatric Pulmonol. 2007;42 :723-728.
- **10- Martinez FD ,Wright AL, Taussing LMHolberg CJ ,Halonen M Morgan WJ ;** Asthma and wheezing in 1st 6 years of life ,The he Group Health Medical Associate N. Eng. J. Med 1995,332:133-138

- 11- James A, Hammar KS, Reinius L, Konradsen JR, Dahl_en S-E, Cilla S€oderh€all C & Hedlin G, A longitudinal assessment of circulating YKL-40 levels in preschool children with wheeze Pediatr Allergy Immunol. 2017;28(1):79-85
- **12- Jang YY & Park JH and Chung HL:** Serum YKL-40 levels may help distinguish exacerbation of post-infectious bronchiolitis obliterans from acute bronchiolitis in young children. Eur J Pediatr (2017) 176:971–978
- **13-Marguet C, Lubrano M, Gueudin M, Le Roux P, Deschildre A, Forget C et al,** In very young infants severity of acute bronchiolitis depends on carried viruses. PLoS One 4 (2):e4596 ,1-6
- **14- Bara I, Ozier A, Girodet PO, et al.**, Role of YKL-40 in bronchial smooth muscle remodeling in asthma. Am J Respir Crit Care Med 2012: 185: 715–22.
- **15 Konradsen JR, James A, Nordlund B, et al.** ,The chitinase-like protein YKL-40; a possible biomarker of inflammation and airway remodeling in severe pediatric asthma. J Allergy Clin Immunol 2013: 132: 328–35.
- **16- Chupp GL, Lee CG, Jarjour N, et al.**, A chitinaselike protein in the lung and circulation of patients with severe asthma. N Engl J Med 2007: 357: 2016– 27.
- **17- Oommen A, Patel R, Browning M, Grigg J.**, Systemic neutrophil activation in acute preschool viral wheeze. Arch Dis Child 2003:88: 529–31.
- **18.** Moldoveanu B, Otmishi P, Jani P, Walker J, Sarmiento X, Guardiola J, et al., Inflammatory mechanisms in the lung. J Inflamm Res 2009; 2:1– 11.
- **19.** Quinton LJ, Jones MR, Robson BE, Simms BT, Whitsett JA, Mizgerd JP., Alveolar epithelial STAT3, IL-6 family cytokines, and host defense during Escherichia coli pneumonia. Am J Respir Cell Mol Biol 2008; 38:699–706.
- **20- Chupp GL, Lee CG, Jarjour N et al (2007),** A chitinase-like protein in the lung and circulation of patients with severe asthma. N Engl J. Med 357:2016–27
- **21. Furuhashi K, Suda T, Nakamura Yet al (2010)** Increased expression of YKL-40, a chitinase-like protein, in serum and lung of patients with idiopathic pulmonary fibrosis. Respir Med 104:1204–10
- **22. Gumus A, Kayhan S, Cinarka H et al (2013),** High serum YKL-40 level in patients with COPD is related to hypoxemia and disease severity. Tohoku J Exp Med 229:163–170

- 23. Nordenbaek C, Johansen JS, Junker P, Borregaard N, Sorensen O, Price PA., YKL-40, a matrix protein of specific granules in neutrophils, is elevated in serum of patients with communityacquired pneumonia requiring hospitalization. J Infect Dis 1999:180: 1722–6.
- **24.** Volck B, Price PA, Johansen JS, et al., YKL- 40, a mammalian member of the chitinase family, is a matrix protein of specific granules in human neutrophils. Proc Assoc Am Physicians 1998: 110: 351–60.
- **25. James AJ, Reinius LE, Verhoek M, et al.**, The Chitinase proteins YKL-40 and chitotriosidase are increased in both asthma and COPD. Am J Respir Crit Care Med 2016: 193: 131–42.
- **26. Lee JH, Park KH, Park JW, Hong CS.,** YKL-40 in induced sputum after allergen bronchial provocation in atopic asthma. J Investig Allergol Clin Immunol 2012: 22: 501–7.
- **27. Gavala ML, Kelly EA, Esnault S, et al.**, Segmental allergen challenge enhances chitinase activity and levels of CCL18 in mild atopic asthma. Clin Exp Allergy 2013:43: 187–97.
- 28- El-fouly MA, Bahaa Eldin MA and Abdelaleem Fouly EK, Evaluation of YKL-40 Protein Levels in Asthmatic Children MJMR, 2018, 29, 3, 203-206
- **29- Mohammed IM, Diab SM, Soliman DR and Shalaby WA**, Study of serum YKL-40 in children with bronchial asthma Egyptian Pediatric Association Gazette 2016, 64, 26–31
- **30- Ober, C., Tan Z., Sun Y, Possick JD, Pan L, Nicolae R et al,** Effect of Variation in CHI3L1 on Serum YKL-40 Level, Risk of Asthma, and Lung Function 2008, N Engl J Med 358;16, 1682-91
- **31- James AJ, Reinius LE, i Verhoek M, Gomes A, Kupczyk M, Hammar U**: Increased YKL-40 and Chitotriosidase in Asthma and Chronic Obstructive Pulmonary Disease, American Journal of Respiratory and Critical Care Medicine 2016, 193, 2, 131-142
- **32- Bendary MA, Ibrahim AA, Mustafa AL and El-Shimi OS,** Role of YKL-40 In Patients with Atopic Dermatitis, Benha Journal of Applied Sciences (BJAS) 2018 ,3, 125-126.
- **33- Jeon SH, Lee YS, Yeo IJ, Lee HP, Yoon J, Son DJ, et al,** Inhibition of Chitinase-3-like-1 by K284-6111 Reduces Atopic Skin Inflammation via Repressing Lactoferrin . Immune Netw. 2021 Jun;21(3)1-17

How to cite

ELSAYED, W., Ibrahim, D. Serum Chitinase3-like protein1 (YKL-40) in Wheezy Infants. Zagazig University Medical Journal, 2023; (620-625): -. doi: 10.21608/zumj.2022.161051.2634